

Arkeo Multichannel API

v0.4

21-07-2025

Contents

| | |
|--|----|
| Introduction | 2 |
| Summary of available commands | 2 |
| List of commands | 3 |
| Command: SetActiveChannel | 3 |
| Command: GetActiveChannel | 3 |
| Command: SetChannelSettings | 3 |
| Command: GetChannelSettings | 3 |
| Command: StartChannel | 3 |
| Command: StopChannel | 4 |
| Command: ForceJV | 4 |
| Command: GetChannelState | 4 |
| Command: GetLatestJV | 4 |
| Command: GetIV | 5 |
| Command: GetSensors | 5 |
| The JSON settings string | 6 |
| Example | 6 |
| Description | 7 |
| Notes for settings JSON | 8 |
| The JSON state string | 9 |
| Example | 9 |
| Description | 9 |
| Examples | 10 |
| Python | 10 |
| Changelog | 11 |

Introduction

This API facilitates remote control of the Arkeo Multichannel using TCP on port 6340, enabling users to perform a variety of commands manage measurements.

A command must sent as a JSON string in the form

```
{'command': command, 'parameter': parameter}
```

Each command must be preceded by a 4-byte integer indicating the total length of the command. After each command, a reply is always sent back. Its format depends on the command (see list below). If an unknown command is sent, the string “Not a valid command” is send instead.

Note: Currently only a single client is allowed to connect to the multichannel software.

Summary of available commands

| Command | Description |
|--------------------|---|
| SetActiveChannel | Set the channel upon which all subsequent actions are performed |
| GetActiveChannel | Get the active channel |
| SetChannelSettings | Modify the settings |
| GetChannelSettings | Retrieve the settings |
| StartChannel | Start measurement |
| StopChannel | Stop the measurement process |
| ForceJV | Force a JV measurement |
| GetChannelState | Retrieve the current state |
| GetLatestJV | Get the voltage and current of last performed JV |
| GetIV | Get the live voltage and current of all channels |
| GetSensors | Get the live sensor values of all connected sensors |

List of commands

Command: SetActiveChannel

Description: Set the channel upon which all succeeding actions are performed.

Parameters:

- channel_id (integer): The ID of the channel.

Response:

- channel_id (integer): The ID of the channel.

Command: GetActiveChannel

Description: Get the active channel.

Parameters:

Response:

- channel_id (integer): The ID of the channel.

Command: SetChannelSettings

Description: Set the JSON configuration string for active channel.

Parameters:

- settings (JSON string): The configuration settings in JSON format (see chapter: The JSON settings string).

Response:

- status (string): OK

Command: GetChannelSettings

Description: Retrieve the JSON configuration string of active channel.

Parameters:

Response:

- settings (JSON string): The configuration settings in JSON format (see chapter: The JSON settings string).

Command: StartChannel

Description: Start the measurement process for active channel.

Parameters:

Response:

- status (string): OK

Command: StopChannel

Description: Stop the measurement process for active channel.

Parameters:

Response:

- status (string): OK

Command: ForceJV

Description: Force a JV measurement on active channel if channel is in tracking mode.

Parameters:

- channel_id (integer): The ID of the active channel.

Response:

- status (string): OK

Command: GetChannelState

Description: Retrieve the current state of a specific channel (e.g., running, stopped, JV, tracking).

Parameters:

Response:

- state (string): JSON of the state of the active channel (see chapter: [The JSON state string](#))

Command: GetLatestJV

Description: Retrieve the voltage and current values of the latest JV.

Parameters:

Response:

- JV Data (string):
v_fw1|j_fw 1|...|v_fw n|j_fw n|| v_rv1|j_rv 1|...|v_rv n|j_rv n

Note:

Forward and Reverse scans are separated by || (double pipe character)
Voltage and Current arrays are interleaved and separated by | (single pipe character)
Voltage unit: V

Current unit A/cm²

Command: GetIV

Description: Retrieve the live voltage and current values of all channels.

Parameters:

Response:

- Voltage and Current values (string):
v1|j1|...|vn|jn

Note:

Voltage and Current arrays are interleaved and separated by | (single pipe character)

Voltage unit: V

Current unit A/cm²

Command: GetSensors

Description: Retrieve the live voltages of all sensors.

Parameters:

Response:

- Sensor values (string):
v1|...|vn|

Note:

Sensor voltages are separated by | (single pipe character)

Voltage unit: V

The JSON settings string

The settings of a channel are represented in a JSON string. When using the **Command: SetChannelSettings**, this is value that the software expects. It is recommended to use the **Command: GetChannelSettings** to read the actual settings and modify those settings.

Example

Below is an example of a channel settings string. A detailed description can be found on the next page.

```
{
  "Index": "1A",
  "Enable": true,
  "User": "Cicci Research",
  "Device": "Si cell",
  "Channel": {
    "VoltageLimit": "10 V",
    "CurrentLimit": 0,
    "InvertedStructure": false
  },
  "JV": {
    "Vmin (V)": -0.1,
    "Vmax (V)": 1.2,
    "Step (mV)": 20,
    "ScanRate (mV/s)": 100,
    "VocDetect": true,
    "Overvoltage (%)": 0,
    "ScanOrder": "FW then RV"
  },
  "Tracking": {
    "TrackEnable": true,
    "Algorithm": "MPPT",
    "Perturbation (V)": 0.02,
    "ConstantOutput": 0,
    "SaveInterval (s)": 10,
    "jvInterval": {"Value": 10, "Unit": "min"},
    "TestDuration": {"Value": 100, "Unit": "hours"}
  },
  "Cell": {
    "Type": "Cell",
    "Area (cm2)": 1,
    "NrCells": 1,
    "NrW cells": 1,
    "W-cellArea (cm2)": 1
  },
  "Note": ""
}
```

Description

| PARAMETER | DESCRIPTION | EXAMPLE | UNIT | TYPE |
|----------------------------------|---|------------|-----------------|---------|
| Index | Unique channel identifier | 1A | | string |
| Enable | Enables or disables the channel | true | | boolean |
| User | Name of the user | User | | string |
| Device | Device Name | Sample | | string |
| Channel | | | | object |
| - VoltageLimit ¹ | Maximum Voltage | 10 V | | enum |
| - CurrentLimit | Current Range | 0 | | integer |
| - InvertedStructure | Inverts the applied voltage | false | | boolean |
| JV | | | | object |
| - Vmin (V) | Minimum voltage | -0.1 | V | float |
| - Vmax (V) | Maximum voltage | 2 | V | float |
| - Step | Voltage increment per step | 20 | mV | integer |
| - ScanRate (mV/s) | Rate at which the voltage is applied | 100 | mV/s | float |
| - VocDetect | Enables the detection of open-circuit voltage | true | | boolean |
| - ScanOrder ² | Order of scanning | FW then RV | | enum |
| Tracking | | | | object |
| - TrackEnable | Enables tracking | true | | boolean |
| - Algorithm ³ | Specifies the algorithm used for tracking | MPPT | | enum |
| - Perturbation (V) | Voltage differential for the tracking algorithm | 0.01 | V | float |
| - ConstantOutput | Setting to maintain a constant output | 0.2 | * | float |
| - SaveInterval (s) | Time between saved data points | 10 | s | integer |
| - jvInterval.Value | Time between JV scans | 10 | | float |
| - jvInterval.Unit ⁴ | Unit for JV interval | min | | enum |
| - TestDuration.Value | Duration of the tracking test | 100 | | float |
| - TestDuration.Unit ⁴ | Unit for Test Duration | 1 | | enum |
| Cell | | | | object |
| - Type ⁵ | Cell Type | Cell | | enum |
| - Area (cm2) | Area of the cell | 1 | cm ² | float |
| - NrCells | Number of cells | 1 | | integer |
| - NrW cells | Number of W cells | 1 | | integer |
| - W-cellArea (cm2) | The area of each W cell | 1 | cm ² | float |
| Note | Free to use field for comments | | | string |

Notes for settings JSON

Detailed explanation and enums of the settings JSON

1. Voltage Limit

The voltage limit of each board is 10 V. 20 V can be reached by connecting 2 boards in parallel. Connect the positive contacts of 2 boards to your device and the negative contacts to each other.

2. Scan Order enum

- 0: FW then RV
- 1: RV then FW
- 2: Forward Only
- 3: Reverse Only

3. Tracking Algorithm enum

- 0: Open circuit
- 1: Short circuit
- 2: MPPT
- 3: MPPT-Stab
- 4: MPPT INC
- 5: Fixed Voltage
- 6: Fixed Voltage (no track)
- 7: Fixed Current
- 8: JV

4. Time unit enum

- 0: seconds
- 1: minutes
- 2: hours

5. Cell Type enum

- 0: Cell
- 1: Parallel Module
- 2: Z Module
- 3: W Module

The JSON state string

The state of a channel is represented in a JSON string. This string is read-only and is obtained using the **Command: GetChannelState**.

Example

Below is an example of a state string. A detailed description can be found below.

```
{ "Enable": false,
  "Channel": "1A",
  "User": "User",
  "Measurement": "JV",
  "Direction": "Forward",
  "State": "Running",
}
```

Description

| PARAMETER | DESCRIPTION | EXAMPLE | DATA TYPE |
|-------------|---------------------------------|---------|-----------|
| Enable | Enables or disables the channel | true | Boolean |
| Channel | ID of the channel | 1A | String |
| User | Name of the user | User | String |
| Measurement | State of the Measurement | JV | String |
| Direction | Current direction of the JV | Forward | String |
| State | Current state of the channel | Running | String |
| | – Idle | | |
| | – Ready to start | | |
| | – Running | | |
| | – Paused | | |
| | – Stopped | | |
| | – Error | | |

Examples

Python

```
import socket
import json

def exampleCommand(command, data=''):
    with socket.socket(socket.AF_INET, socket.SOCK_STREAM) as sock:
        sock.connect(('localhost', 6340))

        message = json.dumps({'command': command, 'data': data})

        # Prepare the length of the string
        length = len(message)
        length_bytes = length.to_bytes(4, byteorder='big') # 4 bytes to
represent the length

        sock.sendall(length_bytes) # Send the length of the string
        sock.sendall(message.encode('utf-8')) # Send the string

        # Receive response
        response = sock.recv(1024)
        print('Received:', response.decode('utf-8'))

exampleCommand('GetChannelSettings')
```

Changelog

v0.1

Initial version

v0.2

Added ForceJV and GetChannelState

v0.3

Added GetLatestJV, GetIV and GetSensors

v0.4

Reorganized the settings JSON and clarified the documentation

Added JSON schema validation